



Nitrite cuvette test 2 - 90 mg/L NO₂-N, 25 tests

Product #:

LCK343

Hazardous

AED Price:

Contact Hach



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	EN ISO 26777, DIN 38405 D10
Instrument:	DR3900, DR6000, DR1900
Measuring range:	2 - 90 mg/L NO ₂ -N
Method:	Diazotization
Number of tests:	25
Parameter:	Nitrite
Platform:	LCK
Storage conditions:	15 - 25 °C