



## Anammox Activity cuvette test, 0 - 1000 mAbs, 25 tests

**Product #:** LCK411.00  
**AED Price:** Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

The test kit LCK411.00 includes one box of LCK411 (pk/25) and one box of LCW906 (pk/25).

LCK411 requires additional accessories that need to be ordered separately: SM10311841 (Filters, 5 - 12 µm particle retention, pk/100) and 2187569 (Porphyrine 2 Reagent, pk/100).

Optional accessories: SM106404 (Funnels, pk/24), SM212864457 (Baffled Flask, pk/1)

### Precise and reliable measurement values

### Maximum safety for users, thanks to the closed system and low amounts of reagents

### Convenient and error-free dosing of the reagents

### Barcode label for automatic recognition in the photometer

---

## Specifications

Digestion required:	Yes
Instrument:	DR3900, DR6000, DR1900
Measuring range:	0 - 1000 mAbs
Method:	Photometric heme detection
Number of tests:	25
Parameter:	Anammox Activity
Platform:	LCK
Storage conditions:	15 - 25 °C
What's included?:	Kit includes:

- Reagents for 25 tests (LCK411)

- 25 cuvettes (LCW906)

Required accessories (need to be ordered separately):

- Filters, 5 - 12 µm particle retention, pk/100 (SM10311841)

## What's included?

Kit includes: - Reagents for 25 tests (LCK411)- 25 cuvettes (LCW906)Required accessories (need to be ordered separately):- Filters, 5 - 12  $\mu$ m particle retention, pk/100 (SM10311841) - Porphyrine 2 Reagent, pk/100 (2187569)

---

## Required Accessories

- Porphyrin 2 Reagent Powder Pillows, 25 mL, pk/100 (Item 2187569)
- Filters, 70 mm, 5 - 12  $\mu$ m particle retention, pk/100 (Item SM10311841)