



Orthophosphate cuvette test 0.01 - 0.5 mg/L $\text{PO}_4\text{-P}$, 20 tests

Product #:

LCK549

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Convenient handling with 50 mm rectangular cuvettes

Highest accuracy for analysis of ultra low range concentrations

Specifications

According to standard:	ISO 6878-1-1986, DIN 38405 D11-4
Instrument:	DR3900, DR6000
Measuring range:	0.01 - 0.5 mg/L $\text{PO}_4\text{-P}$
Method:	Phosphomolybdenum blue
Number of tests:	20
Parameter:	Phosphate, ortho (trace)
Platform:	LCK
Storage conditions:	15 - 25 °C