



QP1680 Consumables Kit, 5000 analysis

Product #: SMKIT501100
AED Price: Contact Hach

Replacement consumables for approximately 5000 injections, depending on sample

Specifications

What's included?:

- 1x Cotton swabs, Cotton/Plastic, 25 pcs
 - 2x Furnace injection port, PTFE
 - 2x Content for halogen scrubber, 1 refill
 - 1x Procat combustion tube refill kit, catalyst and quartz filling
 - 1x Procat combustion tube, complete with catalyst filling
 - 1x Filter, Plastic, 5 µm and disposable, 10 pcs
 - 1x O-ring, Silicone, 22 x 3 mm, for furnace tube, 5 pcs
 - 1x O-ring, PTFE, 6.02 x 2.62 mm, for furnace exit, 5 pcs
 - 2x Stirrer bar, PTFE, 15 pcs
 - 2x Furnace injection needle with sleeve, PTFE sleeve/Metal
 - 10x O-ring Injection and IC port,
 - 1x PTFE top seal, for furnace tube/head
-

What's included?

1x Cotton swabs, Cotton/Plastic, 25 pcs
2x Furnace injection port, PTFE
2x Content for halogen scrubber, 1 refill
1x Procat combustion tube refill kit, catalyst and quartz filling
1x Procat combustion tube, complete with catalyst filling
1x Filter, Plastic, 5 µm and disposable, 10 pcs
1x O-ring, Silicone, 22 x 3 mm, for furnace tube, 5 pcs
1x O-ring, PTFE, 6.02 x 2.62 mm, for furnace exit, 5 pcs
2x Stirrer bar, PTFE, 15 pcs
2x Furnace injection needle with sleeve, PTFE sleeve/Metal
10x O-ring Injection and IC port, 1x PTFE top seal, for furnace tube/head