



## Ammonium cuvette test 10-100 mg/L NH<sub>4</sub>-N, 25 tests

**Product #:**

**LCK503**

**AED Price:**

Contact Hach

**Hazardous**



### **LCK cuvette tests: Minimum amounts of reagents, maximum safety.**

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

#### **Precise and reliable measurement values**

**Maximum safety for users, thanks to the closed system and low amounts of reagents**

**Convenient and error-free dosing of the reagents**

**Barcode label for automatic recognition in the photometer and Truecal calibration**

**Differentiation between tests and measuring ranges by means of colour-coding**

---

### **Specifications**

According to standard:	ISO 7150-1, DIN 38406 E5-1, UNI 11669:2017, ISO 23695
Instrument:	DR3900, DR6000, DR1900
Measuring range:	10 - 100 mg/L NH <sub>4</sub> -N
Measuring range (2):	12.8 - 128 mg/L NH <sub>4</sub>
Method:	Indophenol blue
Number of tests:	25
Parameter:	Ammonium
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)
Truecal:	Yes