



EDTA solution for EZ1016 Total Hardness Analyser

Product #:	APPC1016-04
AED Price:	Contact Hach

Ready-to-use Reagents for EZ Series Online Analysers

Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Specifications

Parameter:	Total Hardness
------------	----------------