



QP1680 High-Temperature TOC Laboratory Analyser, with auto sampler, 65 positions

Product #: LPV448.99.00001

AED Price: Contact Hach

The QP1680 is available as a combined TOC and TN_b analyser or for the individual parameters TOC or TN_b.

The most important features include:

- Direct sample injection eliminates sample contact with valves and the built-in injection syringe, which minimizes the risk of sample carry-over.
- Large diameter sample aspiration tubing can handle particles up to 800 µm, expanding possible applications and reducing clogging.
- Integrated stirrer for each sample position homogenizes particle-containing samples before injection.
- Small footprint with integrated 65-position auto-sampler requires less space in the laboratory (an auto-sampler with 96 positions is also available as an alternative).
- Simple operation, data analysis and system diagnosis thanks to an intuitive software package.

Specifications

Ambient temperature:	-20 - 60 °C
Analysis time:	Approx. 3 minutes
Detection limit:	TC, TOC, NPOC: 50 µg/L
Dimensions:	440 mm x 380 mm x 700 mm (H x W x D)
Gas consumption:	150 mL/min*

*Sample preparation for NPOC determination requires additional 300-350 mL/min.

Gas specifications:	Oxygen or synthetic air: minimum 99.998% (4.8) at 3 - 10 bar
Measurement method:	NDIR (non-dispersive Infrared Detection)
Measuring range:	TC, TIC, NPOC: 0 - 30000 mg/L
Model:	QP1680-TOC
Norms:	TOC / NPOC: ASTM D7573, EN 1484, EPA 415.1, EPA 9060A, ISO 8245, SM 5310B, NEN-ISO 20236
Number of sample positions:	65
Operating conditions:	20 - 30 °C; 20 - 80% relative humidity (non-condensing)
Oxidation method:	Catalytic combustion at 680 °C
Parameter:	TOC
Pollution degree:	2
Power:	Analyser: 750 W max.

	PC: 90 W max.
	Monitor: 100 W max.
Power supply:	Analyser: 100 - 240 VAC, 50/60 Hz, 16 A, with protective grounding
	PC: 100 - 240 VAC, 50/60 Hz, 1.6 A, with protective grounding
	Monitor: 100 - 240 VAC, 50/60 Hz, 1.6 A, with protective grounding
Repeatability:	Up to 10 mg/L TC, TIC, NPOC, TN: < 5%
	> 10 mg/L TC, TIC, NPOC, TN: < 2%
Sample volume:	10 - 1000 µL
Temperature:	Furnace temperature max. 1050 °C (depending on configuration)
Warranty:	24 months
Weight:	Approx. 31 kg
What's included?:	TOC Instrument
	Sample tray with 65 positions for 40 mL vials
	Glass sample vials 40 mL (100 pcs)
	USB 2.0 cable A-B
	Set of tubing PTFE blue-stripe 1/8" OD length 2 m including nut & ferrule
	Reagent container with screwcap for UPW, Acid and Waste
	Power cord
	USB storage device with software, test results and Full User Manual
	Printed User Manual

What's included?

TOC Instrument
Sample tray with 65 positions for 40 mL vials
Glass sample vials 40 mL (100 pcs)
USB 2.0 cable A-B
Set of tubing PTFE blue-stripe 1/8" OD length 2 m including nut & ferrule
Reagent container with screwcap for UPW, Acid and Waste
Power cord
USB storage device with software, test results and Full User Manual
Printed User Manual

Required Accessories

- QP1680 TOC/TNb Starter Package (Item SMKIT500000)