



## **Acid solution for EZ4003 and EZ4004 Total Alkalinity Analyser, 0.025 N HCl (10-250 mg/L), 2 L**

**Product #:** APPC4004-01-EU  
**AED Price:** Contact Hach

### **Ready-to-use Reagents for EZ Series Online Analysers**

#### **Ready-to-use**

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### **Confidence in your results**

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### **Less complexity**

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### **Specifications**

Concentration:	0.025 N HCl
Parameter:	Total Alkalinity
Quantity:	2 L