



ULR Sulfite Fluorescence Test Kit, 6-500 µg/L (ppb), 100 Tests

Product #: 34250000
AED Price: Contact Hach

Sulfite Test Kit for use with the DR1300 FL Fluorometer.

Specifications

Measuring range:	6 - 500 µg/L (ppb) SO ₃
Number of tests:	100
Parameter:	Sulfite
What's included?:	Includes 1 bottle Sulfite Reagent 1, 1 bottle Sulfite Reagent 2, 2x 16 mm Sample Vials, Vial Stand

What's included?

Includes 1 bottle Sulfite Reagent 1, 1 bottle Sulfite Reagent 2, 2x 16 mm Sample Vials, Vial Stand