



ULR Total Chlorine Fluorescence Test Kit, 3-100 µg/L (ppb), 100 Tests

Product #:

34252000

Hazardous

AED Price:

Contact Hach



Total Chlorine Test Kit for use with the DR1300 FL Fluorometer.

Specifications

Measuring range:	3 - 100 µg/L (ppb) Cl ₂ total
Number of tests:	100
Parameter:	Total Chlorine
What's included?:	Includes 1 bottle Total Chlorine Reagent 1, 1 bottle Total Chlorine Reagent 2, 2x 16 mm Sample Vials, Vial Stand

What's included?

Includes 1 bottle Total Chlorine Reagent 1, 1 bottle Total Chlorine Reagent 2, 2x 16 mm Sample Vials, Vial Stand