



## Laton Total Nitrogen cuvette test 100-250 mg/L TN, 25 tests

Product #:

LCK438

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer and Truecal calibration

Differentiation between tests and measuring ranges by means of colour-coding

---

### Specifications

According to standard:	EN ISO 11905-1, ISO23697-1
Description:	Total Nitrogen bonded
Digestion required:	Yes
Instrument:	DR3900, DR6000, DR1900
Measuring range:	100 - 250 mg/L TN <sub>b</sub>
Method:	Koroleff digestion (Peroxisulphate), and photometric detection with 2.6-Dimethylphenol
Number of tests:	25
Parameter:	Nitrogen total (Laton)
Platform:	LCK
Standard method:	EN ISO 11905-1 Digestion with Peroxisulphate
Storage conditions:	15 - 25 °C
Truecal:	Yes