



## Sedimeter for Nitratax sc and Uvas sc

<b>Product #:</b>	<b>LZX450</b>
<b>AED Price:</b>	Contact Hach

The sedimenter can be used as a bypass solution. The sample can be let through the sedimenter without prefiltering.

---

### Specifications

Instrument:	Suitable for all Nitratax and Uvas sc sensors
Length:	Intallation plate for wall mount: 800 mm x 495 mm
Manual languages:	Non verbal instruction sheet
Material:	PVC
Pressure range:	Not pressurized, open drain
Weight:	Sedimeter incl. wall plate only: 14.7 kg