



Flow through unit 50 mm for Uvas plus sc probes

Product #: LZX868
AED Price: Contact Hach

Suitable for media with very high TS contents, sewage treatment plant inlet, etc.

Specifications

Cable length:	10 m
Dimensions (H x W x D):	210 mm x 500 mm x 160 mm
Flow rate:	0.5 L/h
Length:	33 mm
Operating temperature range:	2 - 40 °C (sample temperature)
Path Length:	2 mm
Pressure range:	Not pressurized, open drain
Sample input:	Tube ID 4 mm / AD 6 mm
Sample temperature:	2 - 40°C (35.6 - 104°F)
Sample volume:	0.5 - 10 L
Weight:	3.6 kg