



## Singlet Single-Use pH Buffer Kit, pH 7.00 & 10.01, pk/2x10

**Product #:** 2769820  
**AED Price:** Contact Hach  
**Available**

### Singlet Solution Packs

pH Buffer Standard Solution for pH electrode calibration. Sealed foil pouches of pH buffers. pH is  $\pm 0.02$  at 25°C. Pack of 10 of pH 7.00 and 10 of pH 10.01, 20 mL each

#### Ideal for on-site analyses

Sealed pouches are easy to transport - no risk of spillage.

#### Easy to Use

Just open the pouch and insert the electrode.

#### Fresh and reliable

Foil pouches prevent contamination and air exposure.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

---

### Specifications

Color coded?:	Yes
Description:	pH Buffer Calibration reagent
NIST traceable:	Yes
Package Type:	Singlet
pH value(s):	7.00 & 10.01 $\pm 0.02$ pH at 25°C
Quantity:	20 pcs
Special remarks:	Discard solution after use
Storage conditions:	10 - 25 °C (protect from light)