



Buffer Solution, pH 1.09, 500 mL

Product #:

S11M009

Hazardous

AED Price:

Contact Hach

Ships within 2 weeks



pH Buffer Standard Solution for pH electrode calibration with Certificate of Analysis*. Colourless pH buffer, pH is 1.09 ± 0.02 at 25°C solution. 500 mL bottle

*Available for download from uk.hach.com/CoA; using product part number and lot number.

Specifications

Color coded?:	No
Description:	pH Buffer Calibration reagent
NIST traceable:	Yes
Package Type:	Bottle
pH value(s):	1.09 ± 0.02 pH at 25°C
Quantity:	500 mL
Special remarks:	Observed shelf-life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)