



## Buffer Solution, pH 4.65, 500 mL

**Product #:** S11M010  
**AED Price:** Contact Hach  
**Available**

pH Buffer Standard Solution for pH electrode calibration with Certificate of Analysis\*. Colourless pH buffer, pH is  $4.65 \pm 0.02$  at 25°C solution. 500 mL bottle

\*Available for download from [uk.hach.com/CoA](http://uk.hach.com/CoA): using product part number and lot number.

---

### Specifications

|                     |  |
|---------------------|--|
| Color coded?:       | No   |
| Description:        | pH Buffer Calibration reagent  |
| NIST traceable:     | Yes  |
| Package Type:       | Bottle   |
| pH value(s):        | $4.65 \pm 0.02$ pH at 25°C   |
| Quantity:           | 500 mL   |
| Special remarks:    | Observed shelf-life after opening, following recommended storage and use: 3 months |
| Storage conditions: | 10 - 25 °C (protect from light)  |