



Buffer Solution, pH 4.65, 500 mL

Product #: S11M010
AED Price: Contact Hach
Available

pH Buffer Standard Solution for pH electrode calibration with Certificate of Analysis*. Colourless pH buffer, pH is 4.65 ± 0.02 at 25°C solution. 500 mL bottle

*Available for download from uk.hach.com/CoA: using product part number and lot number.

Specifications

Color coded?:	No
Description:	pH Buffer Calibration reagent
NIST traceable:	Yes
Package Type:	Bottle
pH value(s):	4.65 ± 0.02 pH at 25°C
Quantity:	500 mL
Special remarks:	Observed shelf-life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)