



## Total Alkalinity Test Strips, 0-240 mg/L, 50 tests

**Product #:** 2744850  
**AED Price:** Contact Hach  
**Available**

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to drinking water, food, and beverage processors; boiler operators; and aquatic biologists. Contains 50 tests.

**Easy to use, disposable, and inexpensive**

**A great way to obtain quick, quantitative answers in the field or in the lab**

---

### Specifications

Measuring range:	0 - 240 mg/L
Number of tests:	50
Parameter:	Alkalinity - Phenolphthalein
Platform:	Test Strip
Smallest Increments Steps:	Steps: 0, 40, 80, 120, 180, 240 mg/L (ppm)
Test kit model:	Test Strips
Units:	50 tests