



Test strips, Nitrate 0 - 50 mg/L and Nitrite 0 - 3 mg/L, 25 tests

Product #:

2745425

Hazardous

AED Price:
Available

Contact Hach



Ranges

Nitrate NO₃-N: 0 - 50 mg/L

Nitrite NO₂-N: 0 - 3 mg/L

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Specifications

Footnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Measuring range:	Nitrate 0 - 50 mg/L NO ₃ -N mg/L
Measuring range (2):	Nitrite 0 - 3.0 mg/L NO ₂ -N
Parameter:	Nitrate and Nitrite
Platform:	Test Strip
Smallest increment:	Nitrite: 0, 0.15, 0.3, 1, 1.5, 3 mg/L
Smallest Increments Steps:	Nitrate: 0, 1, 2, 5, 10, 20, 50 mg/L
Test kit model:	Test Strips
Units:	25 tests