



## Multi-parameter test strips, 50 pcs

**Product #:**

**2755250**

**Hazardous**

**AED Price:**

Contact Hach

**Available**



5-in-1 test strips for quick, on-site evaluations of free & total chlorine, range 0-10mg/L, in steps 0, 0.5, 1.0, 2.0, 10; total hardness as CaCO<sub>3</sub>, range 0-25mg/L, in steps 0, 1.5, 3, 7, 15, 25; total alkalinity as CaCO<sub>3</sub>, range 0-240mg/L, in steps 0, 40, 80, 120, 180, 240; and pH, range 6.2-8.4, in steps 6.2, 6.8, 7.2, 7.8, 8.4.

**Free Chlorine (0-10 mg/L)**

**Free Chlorine (0-10 mg/L)**

**Total Hardness (as CaCO<sub>3</sub>, 0-25 gpg and 0-425 mg/L)**

**Total Alkalinity (as CaCO<sub>3</sub>, 0-240 mg/L)**

**pH (6.2 - 8.4)**

---

## Specifications

Measuring range:	Various
Parameter:	Multiparameter
Platform:	Test Strip
Units:	50 tests