



Ammonia test strips, 0 - 6.0 mg/L, 25 tests

Product #: 2755325
AED Price: Contact Hach
Available

Measuring range: 0 - 0.25 - 0.50 - 1.0 - 3.0 - 6.0 mg/L; test strips are ideal in pre-test screening to detect the presence of materials that might interfere with lab testing.

Specifications

Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Measuring range:	0 - 6 ppm
Number of tests:	25
Parameter:	Ammonia, Nitrogen, low range - As $\text{NH}_3\text{-N}$, For freshwater
Platform:	Test Strip
Smallest Increments Steps:	Steps: 0, 0.25, 0.5, 1.0, 3.0, 6.0 ppm
Test kit model:	Test Strips
Units:	25 tests