



## Tube, culture, w/ screw cap, 25 x 150 mm, 10/pk

<b>Product #:</b>	<b>2096250</b>
AED Price:	Contact Hach
Available	

The cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

---

### Specifications

Capacity:	Capacity to Neck 50 mL
Diameter:	25 mm
Dimensions (D x L):	25 mm x 150 mm
Length:	150 mm
Quantity:	10 /pk
Special Features:	Tube, 25 x 150 mm