



Calibration Cylinder Kit, Models 1720 E, D and C

Product #:

4415600

Hazardous

AED Price:
Available

Contact Hach



Formazin Calibration Kit for 1720C and 1720D/E Turbidimeter. Includes calibration cylinder, Tensette Pipette, and 500 mL 4000 NTU Formazin.

Specifications

What's included?:

Calibration cylinder, Tensette pipette 1.0-10.0 mL, and 500 mL 4000 NTU Formazin Standard, Instructions.

What's included?

Calibration cylinder, Tensette pipette 1.0-10.0 mL, and 500 mL 4000 NTU Formazin Standard, Instructions.