



Backpack for portable meters, small, nylon textile

Product #:	2984900
AED Price:	Contact Hach

For comfortably transporting portable meters and up to six 50 mL standard solutions.

Convenient

The answer to hands-free portability of analytical testing equipment, these backpacks combine contemporary styling with innovative design.

Customizable

The adjustable divider system allows you to create customized compartments to organize equipment just the way you like it.

Ergonomic

Lightweight, thick padding material gives superior protection during transportation. A stitched foam back panel makes the backpack well suited for either short or long testing adventures. The sternum buckle system adds support when carrying heavier equipment.

A telescoping handle system and a comfort pad provide additional comfort when used as a backpack. The large backpack can also be used as a wheeled bag or a backpack with in-line skate wheels.

Specifications

Dimensions (H x W x D):	381 mm x 305 mm x 152 mm
Height:	ID:
Length:	152 mm
Material:	Nylon, plastic textile
Weight:	Capacity: 9.1 kg