



## ORP buffer solution, 200 mV, 500 mL

|                   |                      |
|-------------------|----------------------|
| <b>Product #:</b> | <b>25M2A1001-115</b> |
| AED Price:        | Contact Hach         |
| Available         |                      |

ORP Reference Solution used for ORP sensor calibration (in resealable plastic bottles).

---

### Specifications

|                  |                          |
|------------------|--------------------------|
| Instrument:      | Digital Differential ORP |
| Measuring range: | 200 mV                   |
| Parameter:       | ORP                      |
| Quantity:        | 500 mL                   |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)