



Nitrate standard, 25 mg/L NO₃⁻ (5.65 mg/L NO₃⁻-N), 500 mL

Product #: LCW828
AED Price: Contact Hach
Ships within 2 weeks

Standard solution for online nitrate determination

Specifications

Concentration:	25 mg/L (5.65 mg/L NO ₃ ⁻ -N)
Instrument:	Nitratax plus sc, Nitratax clear sc, NT3100sc, NT3200sc
Parameter:	Nitrate
Platform:	Solution
Quantity:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)