



## Iron Reagent Powder Pillows, 0.02-3.00mg/L Fe, FerroVer

**Product #:**

**85428**

AED Price:  
Available

Contact Hach

**Hazardous**



Powder reagent in sealed package, for 25mL sample.

FerroVer Iron Reagent: For the determination of total soluble Iron including all ferric ( $\text{Fe}_3^+$ ) and ferrous ( $\text{Fe}_2^+$ ) Iron as well as complexed Iron compounds (Fe-EDTA). Perfect for many applications and a wide range of concentrations. Determines Iron in the presence of Copper.

---

### Specifications

EPA compliant:	Yes
Instrument:	DR/3000, DR/2000
Measuring range:	0.02 - 3.00 mg/L Fe
Method:	FerroVer
Method number:	8008
Number of tests:	1000
Platform:	Powder Pillows
Quantity:	1000 pk

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)