



## FerroVer Iron Reagent Powder Pillows, 0.02-3.00mg/L Fe

**Product #:**

**85499**

AED Price:  
Available

Contact Hach

**Hazardous**



Powder reagent in sealed packages for the determination of Iron by the 1,10 phenantroline method. pk/100. For 25 mL samples.

---

### Specifications

Description:	FerroVer Iron Reagent Powder Pillows, 25mL
Digestion required:	for total iron
EPA compliant:	Yes
Instrument:	DR/3000, DR/2000
Measuring range:	0.02 - 3.00 mg/L Fe
Method:	FerroVer
Method number:	8008
Number of tests:	100
Parameter:	Iron
Platform:	Powder Pillows
Quantity:	100 /pk
Storage conditions:	10 °C - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)