



Iron Reagent Powder Pillows, 0.02-3.00 mg/L Fe

Product #:

2105769

Hazardous

AED Price:
Available

Contact Hach



Powder reagents in sealed packages for determination of Iron, FerroVer method. For 10 mL sample. pk/100.

Specifications

Description:	FerroVer Iron Reagent Powder Pillows, 10mL
Digestion required:	for total iron
EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000
Measuring range:	0.02 - 3.00 mg/L Fe
Method:	FerroVer
Method number:	8008
Number of tests:	100
Parameter:	Iron
Platform:	Powder Pillows
Quantity:	100 /pk
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)