



CuVer 1 Copper Reagent Powder Pillows, 0.04-5.00mg/L Cu, pk/100

Product #:

2105869

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed packages for the determination of Copper, Bicinchoninate Method. For 10mL sample.

Specifications

EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000, DR/850, DR/820
Measuring range:	0.04 - 5.00 mg/L Cu
Method:	Bicinchoninate
Method number:	8506
Number of tests:	100
Parameter:	Copper
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)