



ZincoVer 5 Zinc Reagent Powder Pillows

Product #:

2106669

Hazardous

AED Price:
Available

Contact Hach



Powder reagents in sealed package for the determination of Zinc, Zincon Method. For 20mL sample.

Specifications

EPA compliant:	Yes
Measuring range:	0.01 - 3.00 Zn
Method:	Zincon
Method number:	8009
Number of tests:	100
Parameter:	Zinc
Platform:	Powder Pillows
Quantity:	100/pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)