



CuVer 2 Copper Reagent Powder Pillows

Product #:

2188299

Hazardous

AED Price:
Available

Contact Hach



For the determination of copper, 0.04-5.00mg/L Cu, in the presence of high levels of Aluminum, Calcium, and Iron. For 25 mL sample. Pack of 100.

Specifications

Method:	Bichinchoninate
Number of tests:	100
Parameter:	Copper
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)