



Manganese Reagent Set, PAN, 0.006-0.700mg/L Mn

Product #:

2651700

AED Price:
Available

Contact Hach

Hazardous



Self contained dropping bottles and powder reagents in sealed packages for determination of Manganese by the PAN method.

Set contains: Alkal. Cyanide reagent (2122326) 50mL, Ascorbic acid (1457799), PAN-Indicator solution (2122426) 50 mL

Specifications

Instrument:	All except DR/3000, DR/2000, DR/850, DR/820
Measuring range:	0.006 - 0.700 mg/L Mn
Method:	PAN
Method number:	8149
Number of tests:	50
Parameter:	Manganese
Sample volume:	50 mL
Special remarks:	Set contains: Alkal. Cyanide reagent (2122326) 50mL, Ascorbic acid (1457799), PAN-Indicator solution (2122426) 50 mL