



## Silver nitrate digital titrator cartridge, 1.128 N, 1 piece

**Product #:**

**1439701**

AED Price:  
Available

Contact Hach

**Hazardous**



Ready-to-use titrant solution for the analysis of Chloride (Silver Nitrate Method), 10 - 10000 mg/L  $\text{Cl}^-$ .

Approximately 40-100 tests.

### **Accurate results - when and where you need them**

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### **Fast and Flexible**

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### **Versatile**

More than 50 parameter/range combinations are available

### **Economical**

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Specifications

Measuring range:	100 - 10000 mg/L $\text{Cl}^-$
Method:	Silver Nitrate
Method number:	8207
Number of tests:	40 - 100
Parameter:	Chloride
Platform:	Digital Titrator
Titrant:	1.128 N $\text{AgNO}_3$

---

## Required Accessories

- Protective gloves size L, blue, nitrile, powder-free, 50 pieces (Item SM743L)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)