



Cyanuric acid 2 reagent powder pillows, 5-50mg/L, pk/50

Product #: 246066
AED Price: Contact Hach
Available

Powder reagent in sealed packages for determination of cyanuric acid by the melamine turbidimetric method. For 25mL sample.

Specifications

Instrument:	See Methods/Procedures for instruments compatible with this product
Measuring range:	5 - 50 mg/L
Method:	Turbidimetric
Method number:	8139
Number of tests:	50
Parameter:	Cyanuric acid
Platform:	Powder Pillows
Quantity:	50 /pk

Required Accessories

- Protective gloves size L, blue, nitrile, powder-free, 50 pieces (Item SM743L)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)