



Cyanuric acid 2 reagent powder pillows, 5-50mg/L, pk/50

Product #: 246066
AED Price: Contact Hach
Available

Powder reagent in sealed packages for determination of cyanuric acid by the melamine turbidimetric method. For 25mL sample.

Specifications

Instrument:	All except DR 6000, DR/4000, PC II
Measuring range:	5 - 50 mg/L
Method:	Turbidimetric
Method number:	8139
Number of tests:	50
Parameter:	Cyanuric acid
Platform:	Powder Pillows
Quantity:	50 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)