



Digital titrator cartridge, sulfuric acid, 1.600 N, 1 piece

Product #:

1438901

Hazardous

AED Price:
Available

Contact Hach



Ready-to-use titrant solution for acid-base and alkalinity titrations.

Specifications

EPA compliant:	Yes
Measuring range:	100 - 4000 mg/L CaCO_3
Measuring range (2):	1 - 10 meq/L Base
Method:	Punto final de fenolftaleína
Method Number 1 Chemistry:	8203
Number of tests:	40 - 100
Parameter:	Alkalinity, Total
Platform:	Digital Titrator
Titrant:	1.600 N H_2SO_4

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)