



Free Chlorine DPD Reagent Powder Pillows, 10 mL, pk/100

Product #:

2105569

Hazardous

AED Price:
Available

Contact Hach



For free chlorine determination by the DPD Free Chlorine method. USEPA accepted for reporting for drinking water analysis using Hach Method 8021.

Range: 0.02 - 2.00 mg/L Cl_2 . Sample Size: 10 mL.

Pack of 100 powder pillows.

Boric Acid free!

Specifications

Description:	DPD Free Chlorine Reagent Powder Pillows, 10mL
EPA compliant:	Yes
Instrument:	All except DR/3000, DR/2000
Measuring range:	0.02 - 2.00 mg/L Cl_2
Method:	DPD
Method number:	8021
Number of tests:	100
Parameter:	Chlorine, free
Platform:	Powder Pillows
Sample size:	10 mL
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)