



Glycine reagent solution, 10 %, 29 mL SCDB

Product #:	2762133
AED Price:	Contact Hach
Available	

For the determination of Chlorine Dioxide, DPD/Glycine Method, 0.04 - 5.00 mg/L ClO_2 .

Glycine Reagent for Hach Chlorine Dioxide DPD methods. 29 mL Self-Contained Dropping Bottle

Specifications

Method:	DPD
Method number:	8065
Number of tests:	100
Platform:	Solution
Quantity:	29 ml

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)