



DPD Total Chlorine Reagent Powder Pillows, 5 mL, pk/100

Product #:

1407699

Hazardous

AED Price:
Available

Contact Hach



For determination of total chlorine by the DPD method. For use with Hach's Chlorine Visual Test Kits. Sample Size: 5 mL. Range: up to 3.5 mg/L Cl_2 . Pack of 100 powder pillows.

Boric Acid free!

Specifications

Method:	DPD
Number of tests:	100
Parameter:	Chlorine, Bromine, Iodine
Platform:	Powder Pillows
Sample size:	5 mL
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)