



EDTA Digital Titrator Cartridge, 0.714 M

Product #:

1495901

Hazardous

AED Price:
Available

Contact Hach



Ready-to-use titrant solution for the analysis of Total Hardness, 10-100+ G.d.H. CaCO_3 .

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Measuring range:	10 - 200 G.d.h.
Method:	EDTA (CalVer® 2)
Method number:	8204
Number of tests:	40 - 100
Parameter:	Hardness, Calcium
Platform:	Digital Titrator
Titrant:	0.714 M EDTA

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)