



## Chloride Reagent Set, 0.1-25.0 mg/L $\text{Cl}^-$

**Product #:**

**2319800**

**Hazardous**

**AED Price:**

Contact Hach

**Ships within 2 weeks**



Reagent set for the determination of Chloride, Mercuric Thiocyanate method.

---

### Specifications

Footnote:	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Instrument:	All except DR/800, PC II
Measuring range:	0 - 25.0 mg/L $\text{Cl}^-$
Method:	Mercuric Thiocyanate
Method number:	8113
Number of tests:	100
Parameter:	Chloride
Platform:	Solution

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)