



Silica Reagent Set, 0.010-1.600 mg/L SiO₂

Product #:

2346800

Hazardous

AED Price:
Available

Contact Hach



For the determination of Silica, Heteropoly Blue Method. For 25 mL sample. 100 tests

Specifications

Instrument:	DR/2000
Measuring range:	0.010 - 1.600 mg/L SiO ₂
Method:	Heteropoly Blue
Method number:	8186
Number of tests:	100
Parameter:	Silica
Sample size:	25 mL
Special Features:	LR, 25 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)