



Silica Reagent Set, 0.010-1.600 mg/L SiO₂

Product #:

2459300

Hazardous

AED Price:
Available

Contact Hach



For the determination of Silica, Heteropoly Blue Method. For 10 mL sample.

Specifications

Instrument:	All except DR/3000, DR/2000, DR/820, PC II
Measuring range:	0.010 - 1.600 mg/L SiO ₂
Method:	Heteropoly Blue
Method number:	8186
Number of tests:	100
Parameter:	Silica
Sample size:	10 mL
Special Features:	LR, 10 mL, 100 tests

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)