



DEHA Reagent 1 Powder Pillows, 3-450µg/L DEHA

Product #:	2167969
AED Price:	Contact Hach
Available	

Powder reagents in sealed packages for the determinaiton of oxygen scavengers by the iron reduction method. pk/100.

Specifications

Description:	DEHA Reagent 1 Powder Pillows
Measuring range:	3 - 450 µg/L DEHA
Measuring range (2):	5 - 600 µg/L Carbohydrazide
Measuring range (3):	9 - 1000 µg/L Hydroquinone
Method:	Iron Reduction
Method number:	8140
Number of tests:	50
Parameter:	Oxygen Scavengers
Platform:	Powder Pillows
Quantity:	100/pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)