



COD TNT, 20 - 1500 mg/L O₂

Product #: 2125951
AED Price: Contact Hach

Hazardous



Reagent Vials (Ø 16 mm) contain premeasured COD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

Description:	Chemical Oxygen Demand
Instrument:	DR 900, DR 1900
Measuring range:	20 - 1500 mg/L O ₂
Method:	Dichromate
Method number:	8000
Number of tests:	25 (1 cuvette for blank needed)
Parameter:	COD
Platform:	Test 'N Tube
Storage conditions:	15 °C - 25 °C (protect from light)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)