



## BOD Seed Inoculum, Polyseed, pk/50 capsules

**Product #:** 2918700  
**AED Price:** Contact Hach  
**Available**

BOD Seed Inoculum for use in Biochemical Oxygen Demand (BOD5) test procedures. A blend of specialized microbial cultures in an easy-to-use capsule. Designed to provide a uniform standard for the degradation of both industrial and municipal waste in BOD5 analyses. Package of 50 capsules. International Laboratory Supply (InterLab), LTD

- \* Consistent seed source: same bacteria, same type of bacteria, same amount of bacteria blend
- \* Known bacterial content is free from pathogens
- \* Save and convenient pre-measured dry product in easy to use capsules
- \* Two years shelf life

---

### Specifications

Application:	Dilution water preparation
Description:	Seed Source
Formula:	Seed Source
Method:	Dilution
Method number:	8043
Package Type:	Capsules
Parameter:	Biochemical Oxygen Demand
Quantity:	50 pcs

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)