



SulfaVer 4 Sulfate Reagent Powder Pillows, 2-70mg/L SO₄

Product #:

1206599

Hazardous

AED Price:
Available

Contact Hach



Powder reagents in sealed package for determination of sulfate. For 25 mL samples. pk/100

Specifications

Description:	SulfaVer 4 Sulfate Reagent Powder Pillow, 25ml
EPA compliant:	Yes
Footnote:	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Instrument:	DR/4000, DR/3000, DR/2010, DR/2000
Measuring range:	2 - 70 mg/L SO ₄
Method:	SulfaVer 4
Method number:	8051
Number of tests:	100
Parameter:	Sulfate
Platform:	Powder Pillows
Quantity:	100 /pk
Reference method:	USEPA 375.4

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)