



## Iodate-Iodide Digital Titrator Cartridge, 0.3998 N

**Product #:**

**1496101**

**Hazardous**

AED Price:  
Available

Contact Hach



Ready-to-use titrant solution for the analysis of Sulfite, 4 - >400mg/L  $\text{SO}_3$ .

Approximately 40-100 tests.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Specifications

Measuring range:	0 - 800 mg/L $\text{SO}_3^{2-}$
Method:	Iodate-Iodide
Parameter:	Sulfite
Platform:	Digital Titrator
Titant:	0.3998 N $\text{KIO}_3$ -KI

---

## Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)