



Sulfite reagent set 0.1-5.0 mg/L SO_3

Product #: HPT430
AED Price: Contact Hach

Specifications

Measuring range:	0.1 - 5.0 mg/L SO_3
Method:	Hach Method
Number of tests:	100
Parameter:	Sulphite
Special Features:	Sulfite Reagent set
Special remarks:	only for DR 4000, DR 2500, DR 2400 and DR 2010, HACH sample cell (21228), 25x54mm, 10mL, needed
Storage conditions:	15 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)